

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

| | |
|------------------------|--------------------------------|
| Product Code(s) | PL.7076 |
| Product Name | Cryptosporidium Differentiator |
| Pure substance/mixture | Mixture |
| Contains methanol | |

1.2. Relevant identified uses of the substance or mixture and uses advised against

| | |
|----------------------|--------------------------|
| Recommended use | Laboratory chemicals |
| Uses advised against | No information available |

1.3. Details of the supplier of the safety data sheet

| | |
|------------------------|--|
| Importer n/a | Manufacturer Pro-Lab Diagnostics 3 Bassendale Road Bromborough Wirral Merseyside CH62 3QL, U.K. |
|------------------------|--|

For further information, please contact

| | |
|----------------|-------------------------|
| E-mail address | uksupport@pro-lab.co.uk |
|----------------|-------------------------|

1.4. Emergency telephone number

| | |
|---------------------|----------------------|
| Emergency Telephone | +44 (0) 151 353 1613 |
|---------------------|----------------------|

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

| | |
|--|---------------------|
| Acute toxicity - Oral | Category 3 - (H301) |
| Acute toxicity - Dermal | Category 3 - (H311) |
| Acute toxicity - Inhalation (Gases) | Category 3 - (H331) |
| Acute toxicity - Inhalation (Vapours) | Category 3 - (H331) |
| Acute toxicity - Inhalation (Dusts/Mists) | Category 3 - (H331) |
| Specific target organ toxicity — single exposure | Category 1 - (H370) |
| Flammable liquids | Category 2 - (H225) |

2.2. Label elements

Contains methanol



Signal word

Danger

Hazard statements

H301 - Toxic if swallowed
 H311 - Toxic in contact with skin
 H331 - Toxic if inhaled
 H370 - Causes damage to organs
 H225 - Highly flammable liquid and vapour

Precautionary statements

P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking
 P260 - Do not breathe dust/fume/gas/mist/vapours/spray
 P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor
 P321 - Specific treatment (see medical advice on this label)
 P370 + P378 - In case of fire: Use dry chemical, CO₂, water spray or alcohol-resistant foam to extinguish
 P403 + P233 - Store in a well-ventilated place. Keep container tightly closed

Additional information

This product requires child resistant fastenings if supplied to the general public. This product requires tactile warnings if supplied to the general public.

2.3. Other hazards

Causes mild skin irritation.

SECTION 3: Composition/information on ingredients

3.1 Substances

Not applicable

3.2 Mixtures

| Chemical name | Weight-% | EC No (EU Index No) | UK REACH registration number | Classification according to GB CLP (SI 2020/1567 as amended) | Specific concentration limit (SCL) | M-Factor | M-Factor (long-term) |
|--------------------------------|-----------------|---------------------|------------------------------|--|---|----------|----------------------|
| methanol 67-56-1 | 50 - <100% | 200-659-6 | - | Flam. Liq. 2 (H225) Acute Tox. 3 (H301) Acute Tox. 3 (H311) Acute Tox. 3 (H331) STOT SE 1 (H370) | STOT SE 1 :: C>=10% STOT SE 2 :: 3%<=C<10% | - | - |
| hydrochloric acid 7647-01-0 | 0.25 - <0.5% | 231-595-7 | - | Press. Gas () Acute Tox. 3 (H331) Skin Corr. 1A (H314) STOT SE 3 (H335) | Eye Irrit. 2 :: 10%<=C<25% Skin Corr. 1B :: C>=25% Skin Irrit. 2 :: | - | - |

| | | | | | | | |
|--|--|--|--|--|-------------------------------------|--|--|
| | | | | | 10%≤C<25% STOT SE 3 :: C>=10% | | |
|--|--|--|--|--|-------------------------------------|--|--|

Full text of H- and EUH-phrases: see section 16

This product does not contain candidate substances of very high concern at a concentration $\geq 0.1\%$ (UK REACH Article 59)

SECTION 4: First aid measures**4.1. Description of first aid measures**

| | |
|---|--|
| General advice | Show this safety data sheet to the doctor in attendance. Immediate medical attention is required. |
| Inhalation | Remove to fresh air. IF exposed or concerned: Get medical advice/attention. If breathing has stopped, give artificial respiration. Get medical attention immediately. Immediate medical attention is required. Do not use mouth-to-mouth method if victim ingested or inhaled the substance; give artificial respiration with the aid of a pocket mask equipped with a one-way valve or other proper respiratory medical device. If breathing is difficult, (trained personnel should) give oxygen. |
| Eye contact | Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Keep eye wide open while rinsing. Do not rub affected area. Get immediate medical attention. |
| Skin contact | Wash off immediately with soap and plenty of water while removing all contaminated clothes and shoes. Get immediate medical attention. |
| Ingestion | Do NOT induce vomiting. Rinse mouth. Never give anything by mouth to an unconscious person. Get immediate medical attention. |
| Self-protection of the first aider | Remove all sources of ignition. Ensure that medical personnel are aware of the material(s) involved, take precautions to protect themselves and prevent spread of contamination. Use personal protective equipment as required. See section 8 for more information. Do not use mouth-to-mouth method if victim ingested or inhaled the substance; give artificial respiration with the aid of a pocket mask equipped with a one-way valve or other proper respiratory medical device. Do not breathe vapour or mist. |

4.2. Most important symptoms and effects, both acute and delayed

| | |
|-----------------|---|
| Symptoms | Coughing and/ or wheezing. Difficulty in breathing. Prolonged contact may cause redness and irritation. |
|-----------------|---|

4.3. Indication of any immediate medical attention and special treatment needed

| | |
|------------------------|------------------------|
| Note to doctors | Treat symptomatically. |
|------------------------|------------------------|

SECTION 5: Firefighting measures**5.1. Extinguishing media**

| | |
|---------------------------------------|---|
| Suitable Extinguishing Media | Dry chemical. Carbon dioxide (CO ₂). Water spray. Alcohol resistant foam. |
| Large Fire | CAUTION: Use of water spray when fighting fire may be inefficient. |
| Unsuitable extinguishing media | Do not scatter spilled material with high pressure water streams. |

5.2. Special hazards arising from the substance or mixture

Specific hazards arising from the chemical Risk of ignition. Keep product and empty container away from heat and sources of ignition. In the event of fire, cool tanks with water spray. Fire residues and contaminated fire extinguishing water must be disposed of in accordance with local regulations.

5.3. Advice for firefighters

Special protective equipment and precautions for fire-fighters Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear. Use personal protection equipment.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Personal precautions Evacuate personnel to safe areas. Use personal protective equipment as required. See section 8 for more information. Avoid contact with skin, eyes or clothing. Ensure adequate ventilation. Keep people away from and upwind of spill/leak. ELIMINATE all ignition sources (no smoking, flares, sparks or flames in immediate area). Pay attention to flashback. Take precautionary measures against static discharges. All equipment used when handling the product must be grounded. Do not touch or walk through spilled material. Do not breathe vapour or mist.

Other information Ventilate the area. Refer to protective measures listed in Sections 7 and 8.

For emergency responders Use personal protection recommended in Section 8.

6.2. Environmental precautions

Environmental precautions Refer to protective measures listed in Sections 7 and 8. Prevent further leakage or spillage if safe to do so. Prevent product from entering drains.

6.3. Methods and material for containment and cleaning up

Methods for containment Stop leak if you can do it without risk. Do not touch or walk through spilled material. A vapour suppressing foam may be used to reduce vapours. Dyke far ahead of spill to collect run-off water. Keep out of drains, sewers, ditches and waterways. Absorb with earth, sand or other non-combustible material and transfer to containers for later disposal.

Methods for cleaning up Take precautionary measures against static discharges. Dam up. Soak up with inert absorbent material. Pick up and transfer to properly labelled containers.

Prevention of secondary hazards Clean contaminated objects and areas thoroughly observing environmental regulations.

6.4. Reference to other sections

Reference to other sections See section 8 for more information. See section 13 for more information.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Advice on safe handling Use personal protection equipment. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Use grounding and bonding connection when transferring this material to prevent static discharge, fire or explosion. Use spark-proof tools and explosion-proof equipment. Keep in an area equipped with sprinklers. Use according to package label instructions. Handle in accordance with good industrial hygiene and safety practice. Avoid contact with skin, eyes or clothing. Take off contaminated clothing and wash it before reuse. Do not breathe vapour or mist. In case of insufficient

ventilation, wear suitable respiratory equipment. Handle product only in closed system or provide appropriate exhaust ventilation. Do not eat, drink or smoke when using this product.

General hygiene considerations

Do not eat, drink or smoke when using this product. Contaminated work clothing should not be allowed out of the workplace. Regular cleaning of equipment, work area and clothing is recommended. Wash hands before breaks and immediately after handling the product. Avoid contact with skin, eyes or clothing. Wear suitable gloves and eye/face protection. Remove and wash contaminated clothing and gloves, including the inside, before re-use. Do not breathe vapour or mist.

7.2. Conditions for safe storage, including any incompatibilities**Storage Conditions**

Keep containers tightly closed in a dry, cool and well-ventilated place. Keep away from heat, sparks, flame and other sources of ignition (i.e., pilot lights, electric motors and static electricity). Keep in properly labelled containers. Do not store near combustible materials. Keep in an area equipped with sprinklers. Store in accordance with the particular national regulations. Store in accordance with local regulations. Keep out of the reach of children. Store locked up.

7.3. Specific end use(s)

Risk Management Methods (RMM) The information required is contained in this Safety Data Sheet.

SECTION 8: Exposure controls/personal protection**8.1. Control parameters****Exposure Limits**

| Chemical name | United Kingdom |
|--------------------------------|---|
| methanol 67-56-1 | TWA: 200 ppm TWA: 266 mg/m ³ STEL: 250 ppm STEL: 333 mg/m ³ Sk* |
| hydrochloric acid 7647-01-0 | TWA: 1 ppm TWA: 2 mg/m ³ STEL: 5 ppm STEL: 8 mg/m ³ |

Biological occupational exposure limits

This product, as supplied, does not contain any hazardous materials with biological limits established by the region specific regulatory bodies.

Derived No Effect Level (DNEL) - Workers

| Chemical name | Oral | Dermal | Inhalation |
|--------------------------------|------|--|--|
| methanol 67-56-1 | | 20 mg/kg bw/day [4] [6] 20 mg/kg bw/day [4] [7] | 130 mg/m ³ [4] [6] 130 mg/m ³ [4] [7] 130 mg/m ³ [5] [6] 130 mg/m ³ [5] [7] |
| hydrochloric acid 7647-01-0 | | | 8 mg/m ³ [5] [6] 15 mg/m ³ [5] [7] |

Notes

[4]

Systemic health effects.

- [5] Local health effects.
- [6] Long term.
- [7] Short term.

Derived No Effect Level (DNEL) - General Public

| Chemical name | Oral | Dermal | Inhalation |
|--------------------------------|--|--|--|
| methanol 67-56-1 | 4 mg/kg bw/day [4] [6] 4 mg/kg bw/day [4] [7] | 4 mg/kg bw/day [4] [6] 4 mg/kg bw/day [4] [7] | 26 mg/m ³ [4] [6] 26 mg/m ³ [4] [7] 26 mg/m ³ [5] [6] 26 mg/m ³ [5] [7] |
| hydrochloric acid 7647-01-0 | | | 8 mg/m ³ [5] [6] 15 mg/m ³ [5] [7] |

Notes

- [4] Systemic health effects.
- [5] Local health effects.
- [6] Long term.
- [7] Short term.

Predicted No Effect Concentration (PNEC)

| Chemical name | Freshwater | Freshwater (intermittent release) | Marine water | Marine water (intermittent release) | Air |
|---------------------|------------|--------------------------------------|--------------|--|-----|
| methanol 67-56-1 | 20.8 mg/L | 1540 mg/L | 2.08 mg/L | | |

| Chemical name | Freshwater sediment | Marine sediment | Sewage treatment | Soil | Food chain |
|---------------------|-------------------------|--------------------------|------------------|-------------------|------------|
| methanol 67-56-1 | 77 mg/kg sediment dw | 7.7 mg/kg sediment dw | 100 mg/L | 100 mg/kg soil dw | |

8.2. Exposure controls

Engineering controls No information available.

Personal protective equipment

Eye/face protection Tight sealing safety goggles.

Hand protection Wear suitable gloves. Impervious gloves.

Skin and body protection Wear suitable protective clothing. Long sleeved clothing. Chemical resistant apron. Antistatic boots.

Respiratory protection No protective equipment is needed under normal use conditions. If exposure limits are exceeded or irritation is experienced, ventilation and evacuation may be required.

General hygiene considerations Do not eat, drink or smoke when using this product. Contaminated work clothing should not be allowed out of the workplace. Regular cleaning of equipment, work area and clothing is recommended. Wash hands before breaks and immediately after handling the product.

Avoid contact with skin, eyes or clothing. Wear suitable gloves and eye/face protection.
Remove and wash contaminated clothing and gloves, including the inside, before re-use. Do not breathe vapour or mist.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

| | |
|-----------------|--------------------------|
| Physical state | Liquid |
| Appearance | Liquid |
| Colour | clear |
| Odour | Alcoholic. |
| Odour threshold | No information available |

| <u>Property</u> | <u>Values</u> | <u>Remarks</u> |
|---|--------------------------|-------------------|
| Melting point / freezing point | ~ -98 °C | ~ -98°C |
| Initial boiling point and boiling range | ~ 65 °C | ~ 65°C @ 1013 hPa |
| Flammability | No data available | None known |
| Flammability Limit in Air | | None known |
| Upper flammability or explosive limits | No data available | |
| Lower flammability or explosive limits | No data available | |
| Flash point | ~ 9.72 °C | None known |
| Autoignition temperature | ~ 455 °C | ~ 455°C |
| Decomposition temperature | | None known |
| pH | No data available | None known |
| pH (as aqueous solution) | No data available | None known |
| Kinematic viscosity | No data available | None known |
| Dynamic viscosity | No data available | None known |
| Water solubility | Soluble in water | None known |
| Solubility(ies) | No data available | None known |
| Partition coefficient | log Pow: ~ -0.77 | None known |
| Vapour pressure | ~ 130.3 hPa @ 2°C | None known |
| Relative density | ~ 0.8 | None known |
| Bulk density | No data available | |
| Liquid Density | No data available | |
| Relative vapour density | ~ 1.11 | None known |
| Particle characteristics | | |
| Particle Size | No information available | |
| Particle Size Distribution | No information available | |
| Explosive properties | No information available | |
| Oxidising properties | No information available | |

9.2. Other information

SECTION 10: Stability and reactivity

10.1. Reactivity

Reactivity No information available.

10.2. Chemical stability

Stability Stable under normal conditions.

Explosion data

Sensitivity to mechanical impact None.

Sensitivity to static discharge Yes.

10.3. Possibility of hazardous reactions

Possibility of hazardous reactions None under normal processing.

10.4. Conditions to avoid

Conditions to avoid Heat, flames and sparks. Excessive heat.

10.5. Incompatible materials

Incompatible materials None known based on information supplied.

10.6. Hazardous decomposition products

Hazardous decomposition products None known based on information supplied.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Information on likely routes of exposure

Product Information

| | |
|---------------------|---|
| Inhalation | Specific test data for the substance or mixture is not available. Toxic by inhalation. (based on components). |
| Eye contact | Specific test data for the substance or mixture is not available. |
| Skin contact | Specific test data for the substance or mixture is not available. Toxic in contact with skin. (based on components). Causes mild skin irritation. |
| Ingestion | Specific test data for the substance or mixture is not available. Toxic if swallowed. (based on components). |

Symptoms related to the physical, chemical and toxicological characteristics

Symptoms Coughing and/ or wheezing. Difficulty in breathing. Prolonged contact may cause redness and irritation.

Acute toxicity

Numerical measures of toxicity

The following values are calculated based on chapter 3.1 of the GHS document

| | |
|-------------------------------|--------------|
| ATEmix (oral) | 101.00 mg/kg |
| ATEmix (dermal) | 303.00 mg/kg |
| ATEmix (inhalation-gas) | 707.0707 ppm |
| ATEmix (inhalation-dust/mist) | 0.504 mg/l |
| ATEmix (inhalation-vapour) | 3.03 mg/l |

Unknown acute toxicity

Component Information

| Chemical name | Oral LD50 | Dermal LD50 | Inhalation LC50 |
|-------------------|-------------------------|--------------------------|-------------------------|
| methanol | = 6200 mg/kg (Rat) | = 15840 mg/kg (Rabbit) | = 22500 ppm (Rat) 8 h |
| hydrochloric acid | 238 - 277 mg/kg (Rat) | > 5010 mg/kg (Rabbit) | = 1.68 mg/L (Rat) 1 h |

| | | | |
|--|--|--|--|
| | | | |
|--|--|--|--|

Delayed and immediate effects as well as chronic effects from short and long-term exposure

| | |
|--|---|
| Skin corrosion/irritation | Classification based on data available for ingredients. Causes mild skin irritation. |
| Serious eye damage/eye irritation | No information available. |
| Respiratory or skin sensitisation | No information available. |
| Germ cell mutagenicity | No information available. |
| Carcinogenicity | No information available. |
| Reproductive toxicity | No information available. |
| STOT - single exposure | Based on the classification criteria of the Globally Harmonized System as adopted in the country or region with which this safety data sheet complies, this product has been determined to cause systemic target organ toxicity from acute exposure. (STOT SE). Causes damage to organs if swallowed. Causes damage to organs in contact with skin. Causes damage to organs if inhaled. |
| STOT - repeated exposure | No information available. |
| Aspiration hazard | No information available. |
| Other adverse effects | No information available. |

SECTION 12: Ecological information

12.1. Toxicity

Ecotoxicity

Unknown aquatic toxicity Contains 0 % of components with unknown hazards to the aquatic environment.

| Chemical name | Algae/aquatic plants | Fish | Toxicity to microorganisms | Crustacea |
|---------------|----------------------|---|----------------------------|-----------|
| methanol | - | LC50: =28200mg/L (96h, Pimephales promelas) LC50: >100mg/L (96h, Pimephales promelas) LC50: 19500 - 20700mg/L (96h, Oncorhynchus mykiss) LC50: 18 - 20mL/L (96h, Oncorhynchus mykiss) LC50: 13500 - 17600mg/L | - | - |

| | | | | |
|--|--|----------------------------|--|--|
| | | (96h, Lepomis macrochirus) | | |
|--|--|----------------------------|--|--|

12.2. Persistence and degradability

Persistence and degradability No information available.

12.3. Bioaccumulative potential

Bioaccumulation

Component Information

| Chemical name | Partition coefficient |
|---------------|-----------------------|
| methanol | -0.77 |

12.4. Mobility in soil

Mobility in soil No information available.

12.5. Results of PBT and vPvB assessment

PBT and vPvB assessment No information available.

| Chemical name | PBT and vPvB assessment |
|-------------------|---------------------------------|
| methanol | The substance is not PBT / vPvB |
| hydrochloric acid | The substance is not PBT / vPvB |

12.6. Endocrine disrupting properties

No information available.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste from residues/unused products Should not be released into the environment. Dispose of in accordance with local regulations. Dispose of waste in accordance with environmental legislation.

Contaminated packaging Empty containers pose a potential fire and explosion hazard. Do not cut, puncture or weld containers.

SECTION 14: Transport information

IATA

- 14.1 UN number or ID number UN1230
- 14.2 UN proper shipping name Methanol
- 14.3 Transport hazard class(es) 3
- Subsidiary hazard class 6.1
- 14.4 Packing group II
- Description UN1230, Methanol, 3 (6.1), II
- 14.5 Environmental hazards Not applicable
- 14.6 Special precautions for user
- Special Provisions A113
- ERG Code 3L

IMDG

- 14.1 UN number or ID number UN1230

| | |
|--|---|
| 14.2 UN proper shipping name | Methanol |
| 14.3 Transport hazard class(es) | 3 |
| Subsidiary hazard class | 6.1 |
| 14.4 Packing group | II |
| Description | UN1230, Methanol, 3 (6.1), II, (9.7°C c.c.) |
| 14.5 Environmental hazards | Not applicable |
| 14.6 Special precautions for user | |
| Special Provisions | 279 |
| EmS-No | F-E, S-D |
| 14.7 Maritime transport in bulk according to IMO instruments | No information available |

RID

| | |
|-----------------------------------|-------------------------------|
| 14.1 UN number or ID number | UN1230 |
| 14.2 UN proper shipping name | Methanol |
| 14.3 Transport hazard class(es) | 3 |
| Subsidiary hazard class | 6.1 |
| 14.4 Packing group | II |
| Description | UN1230, Methanol, 3 (6.1), II |
| 14.5 Environmental hazards | Not applicable |
| 14.6 Special precautions for user | |
| Special Provisions | 279 |
| Classification code | FT1 |

ADR

| | |
|-----------------------------------|--------------------------------------|
| 14.1 UN number or ID number | UN1230 |
| 14.2 UN proper shipping name | Methanol |
| 14.3 Transport hazard class(es) | 3 |
| Subsidiary hazard class | 6.1 |
| 14.4 Packing group | II |
| Description | UN1230, Methanol, 3 (6.1), II, (D/E) |
| 14.5 Environmental hazards | Not applicable |
| 14.6 Special precautions for user | |
| Special Provisions | 279 |
| Classification code | FT1 |
| Tunnel restriction code | (D/E) |

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulations****Authorisations and/or restrictions on use:**

This product contains one or more substances subject to restriction (UK REACH - Annex XVII).

| Chemical name | Restricted substance per REACH Annex XVII | Substance subject to authorisation per REACH Annex XIV |
|--------------------|---|--|
| methanol - 67-56-1 | Use restricted. See item 69. | - |

Persistent Organic Pollutants

Not applicable

Export Notification requirements

Not applicable

Dangerous substance category per COMAH Regulations 2015 (as amended)

H2 - ACUTE TOXIC

H3 - STOT SPECIFIC TARGET ORGAN TOXICITY - SINGLE EXPOSURE

P5a - FLAMMABLE LIQUIDS

P5b - FLAMMABLE LIQUIDS

P5c - FLAMMABLE LIQUIDS

Named dangerous substances per COMAH Regulations 2015 (as amended)

| Chemical name | Lower-tier requirements (tons) | Upper-tier requirements (tons) |
|-------------------------------|--------------------------------|--------------------------------|
| methanol - 67-56-1 | 500 | 5000 |
| hydrochloric acid - 7647-01-0 | 25 | 250 |

The Ozone-Depleting Substances Regulations 2015

Not applicable

The Biocidal Products Regulations 2001 (as amended)

| Chemical name | The Biocidal Products Regulations 2001 (as amended) |
|-------------------------------|---|
| hydrochloric acid - 7647-01-0 | PT2 |

The Water Environment (Water Framework Directive) (England and Wales) Regulations 2017 (as amended)

Not applicable

Poisons Act 1972 (Explosive Precursors) Regulations (as Amended)

| Chemical name | Poisons and Explosive Precursors |
|-------------------|----------------------------------|
| hydrochloric acid | Poison, Reportable 10 % w/w |

International Inventories**TSCA**

Contact supplier for inventory compliance status

DSL/NDSL

Contact supplier for inventory compliance status

EINECS/ELINCS

Contact supplier for inventory compliance status

ENCS

Contact supplier for inventory compliance status

IECSC

Contact supplier for inventory compliance status

KECL

Contact supplier for inventory compliance status

PICCS

Contact supplier for inventory compliance status

AIIC

Contact supplier for inventory compliance status

NZIoC

Contact supplier for inventory compliance status

Legend:**TSCA** - United States Toxic Substances Control Act Section 8(b) Inventory**DSL/NDSL** - Canadian Domestic Substances List/Non-Domestic Substances List**EINECS/ELINCS** - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances**ENCS** - Japan Existing and New Chemical Substances**IECSC** - China Inventory of Existing Chemical Substances**KECL** - Korean Existing and Evaluated Chemical Substances**PICCS** - Philippines Inventory of Chemicals and Chemical Substances**AIIC** - Australian Inventory of Industrial Chemicals**NZIoC** - New Zealand Inventory of Chemicals**15.2. Chemical safety assessment****Chemical Safety Report**

No information available

SECTION 16: Other information**Key or legend to abbreviations and acronyms used in the safety data sheet**

Full text of H-Statements referred to under section 3

H225 - Highly flammable liquid and vapour
 H301 - Toxic if swallowed
 H311 - Toxic in contact with skin
 H314 - Causes severe skin burns and eye damage
 H331 - Toxic if inhaled
 H335 - May cause respiratory irritation
 H370 - Causes damage to organs

Legend

SVHC: Substances of Very High Concern for Authorisation:

Legend Section 8: Exposure controls/personal protection

| | | | |
|---------|-----------------------------|------|----------------------------------|
| TWA | TWA (time-weighted average) | STEL | STEL (Short Term Exposure Limit) |
| Ceiling | Maximum limit value | * | Skin designation |
| + | Sensitisers | | |

Classification procedure

| Classification according to Regulation (EC) No. 1272/2008 [CLP] | Method Used |
|---|--------------------|
| Acute oral toxicity | Calculation method |
| Acute dermal toxicity | Calculation method |
| Acute inhalation toxicity - gas | Calculation method |
| Acute inhalation toxicity - vapour | Calculation method |
| Acute inhalation toxicity - dust/mist | Calculation method |
| Skin corrosion/irritation | Calculation method |
| Serious eye damage/eye irritation | Calculation method |
| Respiratory sensitisation | Calculation method |
| Skin sensitisation | Calculation method |
| Mutagenicity | Calculation method |
| Carcinogenicity | Calculation method |
| Reproductive toxicity | Calculation method |
| STOT - single exposure | Calculation method |
| STOT - repeated exposure | Calculation method |
| Acute aquatic toxicity | Calculation method |
| Chronic aquatic toxicity | Calculation method |
| Aspiration hazard | Calculation method |
| Ozone | Calculation method |

Key literature references and sources for data used to compile the SDS

Agency for Toxic Substances and Disease Registry (ATSDR)
 U.S. Environmental Protection Agency ChemView Database
 European Food Safety Authority (EFSA)
 European Chemicals Agency (ECHA) Committee for Risk Assessment (ECHA_RAC)
 European Chemicals Agency (ECHA) (ECHA_API)
 EPA (Environmental Protection Agency)
 Acute Exposure Guideline Level(s) (AEGL(s))
 U.S. Environmental Protection Agency Federal Insecticide, Fungicide, and Rodenticide Act
 U.S. Environmental Protection Agency High Production Volume Chemicals
 Food Research Journal
 Hazardous Substance Database
 International Uniform Chemical Information Database (IUCLID)
 National Institute of Technology and Evaluation (NITE)
 Australian National Industrial Chemicals Notification and Assessment Scheme (NICNAS)
 NIOSH (National Institute for Occupational Safety and Health)
 National Library of Medicine's ChemID Plus (NLM CIP)
 National Library of Medicine's PubMed database (NLM PUBMED)
 National Toxicology Program (NTP)
 New Zealand's Chemical Classification and Information Database (CCID)
 Organisation for Economic Co-operation and Development Environment, Health, and Safety Publications
 Organisation for Economic Co-operation and Development High Production Volume Chemicals Programme

Organisation for Economic Co-operation and Development Screening Information Data Set
World Health Organization

Revision date 10/10/2023

**This material safety data sheet complies with the requirements of UK REACH Regulations (SI 2019/758 as amended)
Take note of Directive 98/24/EC on the protection of the health and safety of workers from the risks related to chemical agents at work**

Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

End of Safety Data Sheet